

Leading Thermal Analysis.

**NETZSCH**

## Characterization of Pharmaceuticals and Food by Thermal Analysis Instruments coupled with State-of-the-art Gas Analysis Systems



Dr. J. Blumm, S. Eichholz  
NETZSCH-Gerätebau GmbH  
Wittelsbacherstr. 42  
95100 Selb


- Short Introduction into Thermal Analysis Techniques
- Thermogravimetry (TG) – Principle and Information Output
- Coupling of an Evolved Gas Analyzer to a Thermobalance
  - TG-FTIR Coupling - Impurities in a Pharmaceutical Substance
  - TG-QMS Coupling - Sewage Sludge Contaminations in Soil
  - TG-GC-MS Coupling – Contents of Tobacco Smoke
- Summary

- **Definition (ICTAC):**

Thermal analysis (TA) is a group of techniques in which changes of **physical or chemical properties** of the sample are monitored against **time or temperature**, while the temperature of the sample is programmed.

The temperature program may involve heating or cooling at a fixed rate, holding the temperature constant (isothermal), or any sequence of these.

# Main Methods of Thermal Analysis (TA)



# NETZSCH Thermobalances and Simultaneous Thermal Analyzers





STA 449 *F1 Jupiter*<sup>®</sup>



TG 209 *F1 Libra*<sup>®</sup> - ASC

# TG 209 *F1 Iris*<sup>®</sup> - Technical Design



## **Thermophysical Properties:**

**Density change**

## **Product Identification and Characterization:**

**Mass change**

**Decomposition temperatures**

**Thermal stability**

**Compositional analysis**

**Oxidation behavior**

## **Advanced Material Analysis:**


**Decomposition kinetics**

**Rate-controlled mass loss**


# Analysis of different Carbon Modifications





Fine grained carbon modifications are used in many kinds of different industrial applications. Some of them (soot) are known to be carcinogenic.




# SEM and TEM Measurements on different Carbon Modifications




Carbon black



Ethene soot




Graphite



Diesel soot



# TGA Test Results on on different Carbon Modifications



# Evolved Gas Analysis (EGA) in Thermal Analysis **NETZSCH**

In Thermogravimetry (TG) or Simultaneous Thermal Analysis (STA), the mass change (and transformation energetics) of a sample versus temperature or time is measured.

Evolved Gas Analysis yields additional information regarding the nature (composition) of the gases evolved during a mass-loss step.

In most cases, a **Quadrupole Mass Spectrometer (QMS)**, an **Fourier Transform Infrared Spectrometer (FTIR)** or a **Gas Chromatograph-Mass Spectrometer (GC-MS)** are coupled to a TG/STA system for evolved gas analysis.

## **Possibilities:**

- TG/STA-FTIR Coupling (Transfer Line)
- TG/STA-MS Coupling (Capillary Coupling)
- TG/STA-GC-MS Coupling (Transfer Line)
- STA-MS Coupling (Skimmer)

# Evolved Gas Analysis (EGA) in Thermal Analysis **NETZSCH**

**STA/TG-FTIR Coupling:** Yields information on the composition (absorption bands) of the evolved gases (bonding conditions).

- ⇒ Easy interpretation (spectra data bases) of organic vapors without fragmentation
- ⇒ Symmetrical molecules cannot be detected

**STA/TG-MS Coupling:** Yields information on the composition (mass numbers of elements and molecules) of the evolved gases.

- ⇒ High sensitivity
- ⇒ Fast and easy interpretation of atomic/inorganic vapors and standard gases (H<sub>2</sub>, H<sub>2</sub>O, CO<sub>2</sub>, ....)
- ⇒ Fragmentation, interpretation of organic molecules sometimes difficult


**STA/TG-GC-MS Coupling:** Yields information on the composition (mass numbers of elements and molecules) of the evolved gases.

- ⇒ Very high sensitivity
- ⇒ Separation of the volatiles using the GC column
- ⇒ Interpretation of organic vapors significantly improved
- ⇒ Sometimes slow, special measurement processes need to be used or fast GC systems have to be employed

# TG-FTIR Coupling - Realization




# TG-FTIR Coupling - Realization



# Application: Naftopidil

Naftopidil is used to treat benign prostatic hypertrophy (BPH) which causes urinary problems by prostate gland enlargement. It acts to smooth muscles of the prostate and thus increase urinary flow rate.



**Sample:**

crystalline powder

**Experimental conditions:**

Analytical Instruments:

DSC 204, TG 209-FTIR

Sample mass:

0.5 to 16 mg

Crucible:

Al crucible

Temperature program

up 20 °C...160 °C

Heating rate:


10 K min<sup>-1</sup>

Atmosphere:


N<sub>2</sub>



# Application: Naftopidil



# Application: Naftopidil



# Evolved Gas Analysis (EGA) in Thermal Analysis **NETZSCH**

**STA/TG-FTIR Coupling:** Yields information on the composition (absorption bands) of the evolved gases (bonding conditions).

- ⇒ Easy interpretation (spectra data bases) of organic vapors without fragmentation
- ⇒ Symmetrical molecules cannot be detected

**STA/TG-MS Coupling:** Yields information on the composition (mass numbers of elements and molecules) of the evolved gases.

- ⇒ High sensitivity
- ⇒ Fast and easy interpretation of atomic/inorganic vapors and standard gases (H<sub>2</sub>, H<sub>2</sub>O, CO<sub>2</sub>, ....)
- ⇒ Fragmentation, interpretation of organic molecules sometimes difficult

**STA/TG-GC-MS Coupling:** Yields information on the composition (mass numbers of elements and molecules) of the evolved gases.

- ⇒ Very high sensitivity
- ⇒ Separation of the volatiles using the GC column
- ⇒ Interpretation of organic vapors significantly improved
- ⇒ Sometimes slow, special measurement processes need to be used or fast GC systems have to be employed

STA 449 *F1 Jupiter*<sup>®</sup> - QMS 403  
*Aeolos*<sup>®</sup> Capillary Coupling


**NETZSCH**




# Application: Sewage Sludge Contamination in Soil **NETZSCH**




# TG-MS on Dried Sewage Sludge (I)




# TG-MS on Dried Sewage Sludge (II)



# TG-MS on Contaminated Soil (I)



# TG-MS on Contaminated Soil (II)



# Evolved Gas Analysis (EGA) in Thermal Analysis **NETZSCH**

**STA/TG-FTIR Coupling:** Yields information on the composition (absorption bands) of the evolved gases (bonding conditions).

- ⇒ Easy interpretation (spectra data bases) of organic vapors without fragmentation
- ⇒ Symmetrical molecules cannot be detected

**STA/TG-MS Coupling:** Yields information on the composition (mass numbers of elements and molecules) of the evolved gases.

- ⇒ High sensitivity
- ⇒ Fast and easy interpretation of atomic/inorganic vapors and standard gases (H<sub>2</sub>, H<sub>2</sub>O, CO<sub>2</sub>, ....)
- ⇒ Fragmentation, interpretation of organic molecules sometimes difficult

**STA/TG-GC-MS Coupling:** Yields information on the composition (mass numbers of elements and molecules) of the evolved gases.

- ⇒ Very high sensitivity
- ⇒ Separation of the volatiles using the GC column
- ⇒ Interpretation of organic vapors significantly improved
- ⇒ Sometimes slow, special measurement processes need to be used or fast GC systems have to be employed

Yields information on the composition (mass numbers of elements and molecules) of the evolved gases.

- ⇒ Very high sensitivity
- ⇒ Separation of the volatiles using the GC column
- ⇒ Interpretation of organic vapors can be improved

Anyway GC-MS tests generally require minutes to be carried out. This means that a real time analysis which generates a relationship between time/temperature of the TGA run with the GC-MS is often difficult. Solutions:

- ⇒ **Continuous GC-MS** with limited separation possibilities
- ⇒ **Event controlled GC-MS** triggering (high sensitivity and relation between mass loss step and GC-MS result)


# Fully Integrated TG/STA-GC-MS System

**NETZSCH**



# Coupling with chromatographic pre-separation


**NETZSCH**



# Valve box



Gas sampling



valve off  
 sample gas:  
 1-3-sample box-3-4 exit (pumping)  
 GC purge gas:  
 4-3-UNIS 500 injector-GC column

Gas injection




valve on  
 sample gas:  
 1-3 exit (pumping)  
 GC purge gas:  
 4-3-sample box-3-3-UNIS 500 injector-GC column

# Tobacco – Analysis with TG-GC-MS




- In a coupled TG-GC-MS all atmosphere conditions can be realized and the influence on the evolved gases can be clearly demonstrated
- Small amounts of ground cigarette-tobaccos were analyzed in the TG 209 **F1** *Libra*<sup>®</sup> coupled to the Agilent GC-MS via a heated transfer line, the heated JAS valve box and injector system
- Nitrogen atmosphere was used for the pyrolysis experiments, and Nitrogen with 10.7 % oxygen for the combustion experiments (to provide smooth oxidative conditions); heating rate was always 20 K/min for fast results
- The GC-MS was started by the Proteus software simultaneously with the TG 209 F1 for the continuous experiments, and as well started at the selected DTG threshold values for the event-controlled experiments.

# Cigaret-tobacco: continuous mode TG-DTG-TIC







there is always a precise time-temperature correlation for all results

# Chromatogram evaluation of a continuous run




TIC of injection at 18 min (column temperature 250 °C); the gas separation (retention peaks) is mostly sufficient for the identification of the main fragments and molecules


# Chromatogram evaluation of a continuous run



# Single ion monitoring on time/temperature scale




# Detection of nicotine in tobacco (> 200 °C)



m/z 79 → Pyridine;


m/z 84 → Pyridine, 2-(1-methyl-2-pyrrolidinyl)- [84=main peak of nicotine]

Event-control: tobacco, 1.48 mg, N<sub>2</sub>, 20 K/min



reaching the DTG threshold, the injection of the gas onto the column is started and the TG temperature program goes on hold until the GC-oven program is completed; then the TG heating starts again up to the next threshold

# Event-controlled run: Chromatogram at 300 °C



extract of full chromatogram, 15 peaks identified

# Automatic evaluation of Chromatogram at 300 °C



## NETZSCH Library Search Report

Data File : 034-3-10-21 TOBACCO.D Sample : TOBACCO 1.48 mg

Search Libraries: C:\Database\NIST05.L

Peak #	RT(min)	Area%	Identity	Ref #	CAS #	Hit-Qual
1	1.145	19.86	Heptane, 1-(methylthio)- <i>[background]</i>	21310	020291-61-6	20
2	1.180	6.96	1,3-Pentadiene, (E)-	440	002004-70-8	72
3	1.215	3.93	Furan, 3-methyl-	1145	000930-27-8	72
4	1.232	1.80	3-Vinyl-1-cyclobutene	1066	006555-52-8	49
5	1.272	4.57	Furan, 2,5-dimethyl-	2739	000625-86-5	91
6	1.313	2.63	1,5-Cyclooctadiene, 3,8-dimethyl-	15250	1000155-84-1	64
7	1.347	6.62	1,5-Hexadien-3-yne, 2-methyl-	2421	000820-54-2	74
8	1.423	3.01	alpha.-Bourbonene	59835	1000293-01-9	78
9	1.503	5.63	p-Xylene	4946	000106-42-3	76
10	1.550	2.43	Bicyclo[4.2.0]octa-1,3,5-triene	4758	000694-87-1	81
11	1.613	3.87	2,4-Hexadiene, 3,4-dimethyl-	5877	002417-88-1	38
12	1.723	5.99	Phenol	2533	000108-95-2	90
13	1.775	2.46	Butanoic acid, 3-methyl-, 3,7-dimethyl-2,6-octadienyl ester	83868	1000131-82-9	47
14	1.960	11.44	D-Limonene	15162	005989-27-5	94
15	2.081	3.09	Phenol, 3-methyl-	5233	000108-39-4	94


# Influence of atmosphere

All former test were done in inert atmosphere (nitrogen 50 ml/min); Further tests were done under oxidative atmosphere (nitrogen with 10.7 % oxygen, total 56 ml/min) for demonstration of differences in evolved gases.




Cigarette in 100% Oxygen

# Atmosphere influence – pyrolysis/combustion




The TG and DTG show quite strong changes due to the oxygen influence


# Atmosphere influence – m/z 45 (propylene glycol)



# Atmosphere influence – m/z 96 (2,5-methyl-Furan)



# Atmosphere influence – m/z 108 (4-methyl-Phenol)



- Thermal Analysis Instruments are standard in research and development, e.g. in the organic chemistry.
- Coupling to MS, FTIR or GC-MS yields additional information on the material's behavior at high temperatures and the nature of the evolved gases.
- The integrated TG-GC-MS coupling offers the advantage of separating the evolved gases for easier and more accurate interpretation of organic vapors.
- The different examples clearly demonstrated the benefit of using advanced thermal analysis instruments for the characterization of pharmaceuticals, food and related products.